



1. Before treatment, a wound with a deep cavity after debridement



2. Day 4: The first application of the Kerecis graft is completely absorbed



3. Day 18: The wound responded well to the Kerecis fish skin graft.



5. Day 42: The wound is fully healed.

## Kerecis™ Omega3 Wound

Acellular, reconstructive tissue graft for tissue repair.  
Icelandic, piscine origin, Omega3 rich

### Case Review

Dehisced lower extremity wound  
resolved with  
Kerecis™ Omega3 Wound

#### Patient description

55-year-old male with a lower leg fracture.

Internal fixation on lower leg had to be removed due to an infection.

Despite standard wound treatment after surgery, the wound did not heal for over six months.

#### Treatment

The patient received antibiotics before commencing treatment with Kerecis Omega3.

The medication was continued into the 4th week of treatment with Kerecis Omega3.

Wound healed after 10 weeks of treatment.





Keracis Omega3 is intact fish skin that is rich with naturally occurring Omega3 polyunsaturated fatty acids and is used to regenerate damaged human tissue. When inserted onto damaged human tissue such as a diabetic ulcer, the acellular material recruits the body's cells from the wound perimeter. The cells are then incorporated into the acellular fish skin that then ultimately is converted into functional, living tissue. Studies have shown that cells and stem cells proliferate faster in this structure than in other materials such as mammalian-sourced materials.

Fish skin is similar to human skin in presentation and chemical composition, permeability and electric conductivity. Many patients have described relive from pain when using the product as covering exposed nerve endings can provide pain relief.



## Ordering information

Item no	Description	size	Qty
50200S01B2D	Keracis™ Omega3 Wound	3 x 3.5 cm	10 pcs/box
50200S02B2D	Keracis™ Omega3 Wound	3 x 7 cm	10 pcs/box
50200S03B2D	Keracis™ Omega3 Wound	7 x 10 cm	10 pcs/box

Available in the UK from:  
**Pierson Surgical Ltd**  
 North Bradley House, North Bradley, Trowbridge, BA14 0TA  
[sales@piersonsurgical.com](mailto:sales@piersonsurgical.com)  
 01225 766632  
[www.piersonsurgical.com](http://www.piersonsurgical.com)



**KERACIS**

Keracis LLC  
 8300 Greensboro Drive  
 Suite 800  
 McLean, Virginia 22102  
 United States of America  
 Phone: (703) 287-8752

Keracis Limited  
 Eyrargata 2  
 400 Isafjordur  
 Iceland  
 Phone: +354 562-2601

Email: [info@keracis.com](mailto:info@keracis.com) • [www.keracis.com](http://www.keracis.com)

Complies with United States and  
 European medical device regulations  
 U.S. and international patents granted and pending

